

BLEEDING CONTROL

INSTRUCTIONS

1

GLOVES

FOR YOUR PERSONAL PROTECTION
APPLY THESE FIRST



USE SHEARS TO REMOVE ENOUGH CLOTHING TO EXPOSE WOUND AREA

3

COMPRESSED GAUZE

FOR WOUND PACKING

TO STAUNCH BLEEDING, USE ENOUGH TO FILL WOUND SPACE COMPLETELY. HOLD PRESSURE FOR 3-4 MINUTES.



COVER WITH COMPRESSION WRAP (included 4" Emergency Pressure Dressing).

2

WINDLASS TOURNIQUET

FOR SEVERE BLEEDING

USE ONLY ON ARMS & LEGS. APPLY AT LEAST 2-3 INCHES ABOVE WOUND, NOT ON A JOINT. WRAP AROUND LIMB, SLIDE STRAP THROUGH SLIT AND SECURE TIGHTLY WITH VELCRO. TURN WINDLASS UNTIL BLEEDING STOPS. SECURE WITH CLIP.



4


PRESSURE DRESSING

FOR MODERATE BLEEDING

USE OVER COMPRESSED GAUZE IF WOUND REQUIRES PACKING. PLACE PAD OVER WOUND SITE. WRAP SNUGLY TO APPLY PRESSURE TO WOUND.



SECURE WITH CLIP.



RESCUE ESSENTIALS

www.rescue-essentials.com

312 F Street * PO Box 987 * Salida CO 81201
719-539-4843 info@rescue-essentials.com